

CRS Cytokine Multiplex ELISA Detection Kit

Overview

Cat.No. HG-HC001

The kit is an Enzyme Immunoassay kit for the semi-quantification of Human CAR-T / CRS (Cytokine Release Syndrome) Cytokine (IL2, IL6, IL10, IFN gamma) in serum, plasma and cell culture supernatants.

Specification

Assay range:

- IL2: 15.625-500 pg/ml
- IL6: 31.25-1000 pg/ml
- IL10: 15.625-500 pg/ml
- IFN-γ: 15.625-500 pg/ml

Limit of quantitation:

- IL2: 15.625 pg/ml
- IL6: 31.25 pg/ml
- IL10: 15.625 pg/ml
- IFN-γ: 15.625 pg/ml

Standard curve

